



Serological Pipets

Class A, Reusable Glass, Color-Coded, w/ Colored Markings, To Deliver

SIBATA 2011A-Series Serological Pipets meet ASTM E-1044 Style 1 Specifications.

These pipets are graduated to the tip (tempered) and are also calibrated to deliver. Feature permanent amber markings and color-codings which are enameled onto the glass.

The top end of sizes 5mL to 25mL is constricted.

Specifications

Cat. No.	Capacity (mL)	Tolerance (\pm mL)	Graduation Interval (mL)	Approx. O.D. x Height (mm)	Color Code
2011A-01-c	0.1	0.0025	0.01	6 x 300	White
2011A-02-c	0.2	0.004	0.01	6 x 300	Black
2011A-1-d	1.0	0.01	0.1	7 x 360	Red
2011A-1-c	1.0	0.01	0.01	7 x 360	Yellow
2011A-2-d	2.0	0.01	0.1	7 x 360	Green
2011A-5	5.0	0.02	0.1	8 x 365	Blue
2011A-10	10.0	0.03	0.1	10 x 370	Orange
2011A-25	25.0	0.05	0.1	13 x 455	White

Packaging

Cat. No.	Pack Qty.	Case Qty.
2011A-01-c	12	72
2011A-02-c	12	72
2011A-1-d	12	72
2011A-1-c	12	72
2011A-2-d	12	72
2011A-5	12	72
2011A-10	12	72
2011A-25	12	72



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
 (800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
 www.veegee.com • e-mail: sales@veegee.net