



## Volumetric Pipets

**Class A, Reusable Glass, Color-Coded, w/ Colored Markings, To Deliver**

SIBATA 2040A-Series Volumetric Pipets meet ASTM E-969 Class A Specifications & USP Standards for Volumetric Glassware.

These pipets have sharp and permanent graduation lines and are calibrated to deliver. Feature permanent amber markings and color-codings which are enameled onto the glass.



### Specifications

Cat. No.	Capacity (mL)	Tolerance ( $\pm$ mL)	Approx. Length (mm)	Color Code	Pack Qty.	Case Qty.
2040A-05	0.5	0.005	300	2-Black	12	72
2040A-1	1.0	0.006	345	Blue	12	72
2040A-2	2.0	0.006	355	Orange	12	72
2040A-3	3.0	0.01	360	Black	12	72
2040A-4	4.0	0.01	365	2-Red	12	72
2040A-5	5.0	0.01	375	White	12	72
2040A-6	6.0	0.01	390	2-Orange	12	72
2040A-7	7.0	0.01	410	2-Green	12	72
2040A-8	8.0	0.02	430	Blue	12	72
2040A-9	9.0	0.02	436	Black	12	72
2040A-10	10.0	0.02	440	Red	12	72
2040A-13	13.0	0.03	460	White	6	36
2040A-15	15.0	0.03	460	Green	6	36
2040A-20	20.0	0.03	475	Yellow	6	36
2040A-25	25.0	0.03	480	Blue	6	36
2040A-50	50.0	0.05	540	Red	6	36
2040A-100	100.0	0.08	600	Yellow	6	36



## VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034  
(800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826  
www.veege.com • e-mail: sales@veege.net