



Volumetric Flasks

Class A, Snap Cap, To Contain

SIBATA 2303A-Series Volumetric Flasks meet ASTM E-288 Class A Specifications & USP Standards for Volumetric Glassware.

These flasks feature machine-blown bodies and heavy-beaded, heavy-tubing necks, tooled for snap-cap fit. The graduation line is sharp and permanent; large, white, block-letters allow the labeling to be easily read. Each flask also has a large marking square.


All sizes are supplied with a snap cap.

Specifications

| Cat. No. | Capacity (mL) | Tolerance (\pm mL) | Snap Cap No. | Approx. O.D. x Height (mm) |
|------------|---------------|-----------------------|--------------|----------------------------|
| 2303A-5 | 5 | 0.02 | 1 | 21 x 59 |
| 2303A-10 | 10 | 0.02 | 1 | 26 x 85 |
| 2303A-25 | 25 | 0.03 | 2 | 37 x 110 |
| 2303A-50 | 50 | 0.05 | 3 | 49 x 136 |
| 2303A-100 | 100 | 0.08 | 3 | 58 x 164 |
| 2303A-200 | 200 | 0.10 | 4 | 74 x 205 |
| 2303A-250 | 250 | 0.12 | 4 | 78 x 225 |
| 2303A-500 | 500 | 0.20 | 5 | 96 x 260 |
| 2303A-1000 | 1000 | 0.30 | 6 | 122 x 300 |

Packaging

| Cat. No. | Pack Qty. | Case Qty. |
|------------|-----------|-----------|
| 2303A-5 | 6 | 72 |
| 2303A-10 | 6 | 72 |
| 2303A-25 | 6 | 72 |
| 2303A-50 | 4 | 48 |
| 2303A-100 | 4 | 48 |
| 2303A-200 | 4 | 24 |
| 2303A-250 | 4 | 24 |
| 2303A-500 | 2 | 12 |
| 2303A-1000 | 2 | 8 |



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
 (800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
 www.veegee.com • e-mail: sales@veegee.net