



Micro Volumetric Flasks

Class A, K Stopper, Hexagonal Base, To Contain

SIBATA 2306A-Series Micro Volumetric Flasks meet ASTM E-237 Class A, Style I Specifications.

These flasks feature machine-blown bodies and heavy-beaded, heavy-tubing necks, tooled for K stoppers. The graduation line is sharp and permanent; large, white, block-letters allow the labeling to be easily read.

Both sizes feature hexagonal bases for added stability and are supplied with a glass K stopper.

Specifications

Cat. No.	Capacity (mL)	Tolerance (\pm mL)	K Stopper No.	Approx. O.D. x Height (mm)
2306A-1	1	0.01	9	10 x 79
2306A-2	2	0.015	9	12 x 90

Packaging

Cat. No.	Pack Qty.	Case Qty.
2306A-1	6	96
2306A-2	6	96



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
(800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
www.veegee.com • e-mail: sales@veegee.net