



Glass Stoppers

Standard Taper

SIBATA ∇ Glass Stoppers meet ASTM E-675 Specifications.

These stoppers are lightweight yet very strong to reduce chipping and breakage, thus minimizing replacement costs. The flat top design allows the stopper to stand on its head, which frees the user's hand and minimizes potential contamination.

Fully interchangeable with other common flask stoppers and fit all standard labware of comparable ∇ size.



Specifications

Cat. No.	∇ Stopper No.	Tinting	Approx. Length Of Ground Zone (mm)	Height Above Ground Joint (mm)
03081	9	Clear	14	10
03082	13	Clear	14	13
03083	16	Clear	16.5	18
03084	19	Clear	19	18
03085	22	Clear	20.5	20
03087	27	Clear	21	21
03081-A	9	Amber	14	10
03082-A	13	Amber	14	13
03083-A	16	Amber	16.5	18

Packaging

Cat. No.	Fits 2306A & 2307A Size (mL)	Pack Qty.
03081	5, 10, 25	6
03082	50, 100	6
03083	200, 250	6
03084	500	6
03085	1000	6
03087	2000	6
03081-A	5, 10, 25	6
03082-A	50, 100	6
03083-A	200, 250	6



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
 (800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
 www.veegee.com • e-mail: sales@veegee.net