



## Polyethylene Stoppers

### Standard Taper

SIBATA  $\nabla$  Polyethylene Stoppers meet ASTM E-675 Specifications.

These stoppers are made with closed bottoms and polyethylene, which conforms to  $\nabla$  stopper dimensions. The three unit rings on the base provide an efficient seal.

Can be used in place of conventional  $\nabla$  glass stoppers.

### Specifications

Cat. No.	$\nabla$ Stopper No.
009	9
013	13
016	16
019	19

### Packaging

Cat. No.	Fits 2306A & 2307A Size (mL)	Pack Qty.
009	5, 10, 25	50
013	50, 100	50
016	200, 250	50
019	500	50



## VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034  
(800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826  
www.veege.com • e-mail: sales@veege.net