

# VEE GEE®

## Burets

### Class A, w/ PTFE Stopcocks

VEE GEE 20114-Series Burets meet ASTM E-287 Class A Specifications.

Made from precision, uniform Schott DURAN® true-bore tubing to ensure the stipulated capacity tolerances. These burets feature durable, permanent markings, fine, sharp lines, and large, easy-to-read numbers. Supplied with PTFE stopcocks to eliminate freezing and lubricant contamination and to provide a precision fit within the barrel.

All sizes are supplied with a 2mm stopcock bore plug.

### Specifications

Cat. No.	Capacity (mL)	Graduation Interval (mL)	Tolerance ( $\pm$ mL)	Stopcock Bore (mm)	Approx. O.D. x Height (mm)
20114-10	10	0.05	0.02	2	9 x 510
20114-25	25	0.10	0.03	2	12 x 550
20114-50	50	0.10	0.05	2	13 x 758
20114-100	100	0.20	0.10	2	17 x 775

### Packaging

Cat. No.	Pack Qty.	Case Qty.
20114-10	1	24
20114-25	1	24
20114-50	1	24
20114-100	1	24



## VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034  
(800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826  
www.veegee.com • e-mail: sales@veegee.com