

# VEE GEE®

## Separatory Funnels

Squibb, Pear-Shaped, ⌘ Stopper, Teflon® Stopcock

VEE GEE 20149-Series Separatory Funnels meet ASTM E-1096 Type 4 Specifications for Classification, Design, and Markings.

These separatory funnels feature smooth tooled tips, large, white block letters, and large marking squares. Supplied with Teflon® stopcocks to eliminate freezing and lubricant contamination and to provide a precision fit within the barrel. The stopcock bores, body openings, and stems are carefully aligned for even liquid flow with fast cut-off. Made from Schott DURAN® borosilicate tubing with a low coefficient of expansion. Borosilicate tubing also provides uniform wall thickness resulting in the optimum balance between thermal shock resistance and mechanical strength.

All sizes are supplied with a polyethylene ⌘ stopper.

### Specifications

Cat. No.	Capacity (mL)	⌘ Stopper No.	Stopcock Plug Size	Approx. O.D. x Stem Length (mm)	Approx. Total Length (mm)
20149-250	250	22	2	10 x 65	345
20149-500	500	27	4	10 x 65	392
20149-1000	1000	27	4	12 x 70	400
20149-2000	2000	38	6	12 x 75	570

### Packaging

Cat. No.	Pack Qty.	Case Qty.
20149-250	1	24
20149-500	1	24
20149-1000	1	12
20149-2000	1	9



## VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034  
(800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826  
www.veegee.com • e-mail: sales@veegee.com

Printed in the USA



No. PS112008GL20149BRLA