


BUR

Buret/Retort Stand

Non-corrosive, moulded Retort Stand provides an excellent replacement for metal ones. Even the metallic rod has a polypropylene covering making it non-corrosive.

The rod can be threaded into the heavy base giving way to a perfectly stable retort stand. Interfaces with the VEE GEE Class A Buret Series and Buret Clamp.



Specifications

Cat. No.	Dimensions Base (mm)	Dimensions Rod (mm)	Pack Qty.	Weight
20117	300 X 200	12.5 X 750	1	5.5 LBS. 2.5 KG



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
(800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
www.veegee.com • email: sales@veegee.com

VBRCBuret-Retort121814
Printed in the USA