


PIPETTE

Pipette Pump

Presenting a fast release pipetting device for precise pipetting & emptying. The Knurled thumbwheel can be rotated to draw up or dispense liquids while the releaser may be pushed in for zip-quick emptying.

The moulded chuck has a silicon rubber collet inside to hold various pipettes of standard sizes. It resists acids & alkalis and can be easily disassembled for thorough cleaning.



Specifications

Cat. No.	Capacity (mL)	Pack Qty.
20104-2	2	1
20104-10	10	1
20104-25	25	1



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
(800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
www.veegee.com • email: sales@veegee.com

VBRCPipettePump121914
Printed in the USA